

AGILENT INFINITYLAB SUPPLIES FOR THE AGILENT 1260 INFINITY II LC WITH VIALSAMPLER

Quick reference guide

Agilent supplies for Agilent instruments

Agilent Technologies is committed to optimizing your laboratory's productivity, so we have produced this list of the most commonly ordered supplies and parts for the 1260 Infinity II LC with Vialsampler.

The InfinityLab family is an optimized portfolio of LC instruments, columns, and supplies designed to work together seamlessly for maximum efficiency and performance.



InfinityLab Supplies for Pumps	Part No.
Bottle Head Assembly	G7120-60007
Glass Filter, solvent inlet, 20 µm pore size	5041-2168
Stainless Steel Filter, solvent inlet, 10 µm pore size	01018-60025
Piston Seal PTFE, 2/Pk (reversed phase)	5063-6589
Piston Seal PE, 2/Pk (normal phase)	0905-1420
Piston, sapphire	5063-6586
Passive Inlet Valve	G1312-60066
Passive Inlet Valve (normal phase)	G1312-60166
Active Inlet Valve Body, without cartridge	G1312-60025
Cartridge, for active inlet valve, 600 bar	G1312-60020
Outlet Valve	G1312-60067
Outlet Valve (normal phase)	G1312-60167
Purge Valve	G1312-60071
PTFE Frits, for purge valve (pack of 5)	01018-22707
Seal Cap, for purge valve	5067-4728
Wash Seal PTFE (reversed phase)	0905-1175
Wash Seal PE (normal phase)	0905-1718
Gasket, seal wash (pack of 6)	5062-2484
Peristaltic Pump, for seal wash (silicone tubing)	5067-4793
Peristaltic Pump, for seal wash (pharmed tubing)	5065-4445
Silicone Tubing, seal wash (5 m/pk)	5065-9978

For more information see: www.agilent.com/chem/store

InfinityLab Supplies for Vialsampler	Part No.
Drawer for 66 x 2 mL Vials	G7129-60010
Drawer for 18 x 6 mL Vials	G7129-60110
Drawer-Kit 100 x 2 mL Vials Classic	G7129-68210
Drawer for 50 x 2 mL Vials Classic, left	G7129-60210
Drawer for 50 x 2 mL Vials Classic, right	G7129-60220
Peristaltic Pump, (with Pharmed tubing)	5065-4445
Needle Assembly 1260 Autosampler (600 bar)	G7129-87200
Needle Assembly (slotted) for high injection volumes (400 bar)	G7129-87202
Needle Seat Assembly PEEK 0.17 mm	G7129-87017
Rotor Seal (PEEK) Injection Valve, 600 bar	0101-1416
Sample Loop 100 µL	G7129-60500
Metering Seal 100 µL	0905-1503
Piston Ceramic 100 µL	5067-5678
Sample Loop Extension (900 µL)	G1313-87303
Metering Seal, 900 µL	0905-1294
Piston, sapphire 900 µL	5062-8587

InfinityLab Supplies for Multicolumn Thermostat	Part No.
Quick-Connect Heat Exchanger (ID: 0.17 mm, Vol. 1.6 µL)	G7116-60051
Quick-Connect Heat Exchanger (ID: 0.17 mm, Vol. 3.0 µL)	G7116-60061
Rotor Seal (PEEK) for 4-column switching valve, 600 bar	5068-0264
Rotor Seal (PEEK) for 2 Pos/6 Port valve, 600 bar	0101-1409
Rotor Seal (PEEK) for 2 Pos/10 Port valve, 600 bar	0101-1415
Column ID Tag Assembly	5067-5917
Column Holder Clips, soft (2/pk)	G7116-68003
Column Holder Clamp for larger columns (2/pk)	G7116-68004

InfinityLab Supplies for VWD Detectors	Part No
Deuterium Lamp, with RFID Tag	G1314-60101
Standard Flow Cell 10 mm, 14 µL, with RFID tag	G1314-60186
Semi-micro Flow Cell 6 mm, 5 µL, with RFID tag	G1314-60183
High Pressure Flow Cell 10 mm, 14 µL, 400 bar, with RFID tag	G1314-60182
Inlet Tubing Assembly, PEEK 600 mm x 0.2 mm ID	5062-8522
Waste Tube Kit	5062-8535


InfinityLab Supplies for DAD & MWD G7115A & G7165A	Part No.
Standard Flow Cell, with RFID tag	G1315-60022
Semi-micro Flow Cell, with RFID tag	G1315-60025
Long-life Deuterium Lamp, with RFID tag	2140-0820
Tungsten Lamp	G1103-60001

InfinityLab Supplies for DAD G7117A	Part No.
Long-life Deuterium Lamp (8-pin), with RFID tag	5190-0917
Max-Light Cartridge Cell (10 mm, V(σ) 1.0 µL)	G4212-60008
Max-Light Cartridge Cell (60 mm, V(σ) 4.0 µL)	G4212-60007


InfinityLab Supplies for FLD	Part No.
Flow Cell 8 µL, 20 bar	G1321-60005
Flow Cell 4 µL, 20 bar	G1321-60015
FLD Wavelength Calibration Kit	G7121-68001

InfinityLab Capillaries and Fittings	Part No.
Capillary ST 0.17 mm x 900 mm, Pump to Vialsampler (with ICC)	5500-1217
Capillary ST 0.17 mm x 500 mm, Pump to Sampler & Sampler to MCT	5500-1246
Capillary ST 0.17 mm x 105 mm, Vialsampler to ICC	5500-1240
Capillary ST 0.17 mm x 120 mm, ICC to Column	5500-1250
Capillary ST 0.12 mm x 280 mm, Column to Detector	5500-1191
Quick Connect Fitting Assembly, 0.17 mm ID x 105 mm length	5067-6166
Quick Turn LC Fitting	5067-5966

InfinityLab System Supplies and Kits	Part No.
1260 Infinity II Convenience Kit	5067-6617
1260 Infinity II A-Line Capillary Kit	5067-6614
Inline Filter Kit (0.3 µm)	5067-6189
System Tool Kit	G7120-68708
Solvent Bottle, transparent 1 L	9301-6524
Solvent Bottle, amber 1 L	9301-6526
Solvent Bottle, transparent 1 L, with cap	9301-6528
Solvent Bottle, clear, 500 mL, with cap	9301-6523
Solvent Bottle, amber, 500 mL, with cap	9301-6525
Identification Rings (silicone) for Solvent Bottles (8pcs/4colors)	9301-6529
Sticker for Solvent Bottles (100/pack)	9301-6530
Stay Safe Cap Starter Kit	5043-1222
Stay Safe Cap (GL45) with 1 port, 1x venting valve, 2x fittings 3.2 mm	5043-1217
Venting Valve for Stay Safe Cap	5043-1190
6 L Waste Container including 1x Stay Safe Cap GL45 with 4 Ports	5043-1221
Charcoal Filter (48 g) for Waste Container	5043-1193
Flex Bench	5043-1252
PM Kit Binary pump	G1312-68741
PM Kit Quaternary Pump	G1310-68741
Seal Wash PM Kit for Quaternary and Binary Pumps	G1310-68742



Agilent Technologies
Authorized Distributor

 <http://campilab.by>

000 "KAMNУЛAБ"

For research use only. Not for use in diagnostic procedures.
This information is subject to change without notice.

© Agilent Technologies, Inc. 2016
Printed in USA, May 1, 2016
5991-6888EN